

Immunoglobulin M

- **Immunoglobulin M, total, in serum**

Material: 1 ml serum

Stability: 4 months at 2 to 8°C

TAT: same day

Method: turbidimetry

Units: mg/dl

Ref.- range: see report

- **Immunoglobulin M in CSF (see Reiber Scheme)**

Material: 1 ml CSF

TAT: 7 - 10 days*

Method: turbidimetry

Units: mg/dl

Ref. - range: see report

For complete list of laboratory test offered at Freiburg Medical Laboratory, please visit
<http://www.fml-dubai.com/parameter-listings/>