

# Xylene

## General:

Methylhippuric acid is a metabolite of xylene. Approx. 95 % of the resorbed xylenes are eliminated via the kidneys and detected as methylhippuric acid in urine. The following tests are available:

Material: 10 ml urine

Preamalytics: cooled dispatch

TAT: 7-10 days\*

Method: HPLC

Units: mg/L

Ref.- range: see report

For complete list of laboratory test offered at Freiburg Medical Laboratory, please visit  
<http://www.fml-dubai.com/parameter-listings/>