

Anti staphylolysin antibody

General:

In case of staphylococcus infections (e.g. wound infection), the organism reacts with formation of antibodies against the staphylococcus toxin staphylolysin. An increase of the anti-staphylolysin titer is observed approx. 2-3 weeks after infection.

Indication: Focus search, titer level control

Material: 1 ml serum

TAT: 7-10 days*

Method: PA

Units: IU/ml

Ref.- range: up to 2.0

For complete list of laboratory test offered at Freiburg Medical Laboratory, please visit <http://www.fml-dubai.com/parameter-listings/>