

Endothelial cell autoantibodies

General:

Auto-antibodies against endothelial cells (AECA) include a very heterogeneous auto-antibody group and their specificity is not yet sufficiently clear. Their targets are numerous different endothelial cell antigens, which also occur on other cells and are described with a multiplicity of different clinical symptoms.

Therefore no selective diagnostic relevance can be attributed at the moment to these antibodies. A possible immune pathogenesis of vessel damages might be due to an interaction with endothelial auto-antibodies.

Material: 2 ml serum

TAT: 7-10 days*

Method: IFT

Units: Titer

Ref.- range: <1:100

For complete list of laboratory test offered at Freiburg Medical Laboratory, please visit <http://www.fml-dubai.com/parameter-listings/>