

Urine Culture

Indication: The aim of the test is to isolate, quantify and permit presumptive identification and differentiation of the major microorganisms that could be causing urinary tract infection (UTIs).

Preanalytics: The specimen of urine should be collected aseptically into a sterile universal leak proof container following a disinfection of the perineal area.

Stability: Specimen should be stored and transported in sealed plastic bags. Processing of samples should be done as soon as possible after collection.

Material: 10 ml urine

TAT: 2-3 days, FML

Method: Microscopic /Biochem/Culture /Sensitivity

Units: Qualitative

For complete list of laboratory test offered at Freiburg Medical Laboratory, please visit <http://www.fml-dubai.com/parameter-listings/>