

Methylhippuric acid

General:

Methylhippuric acid is a metabolite of xylene. Approx. 95 % of the resorbed xylenes are eliminated via the kidneys and detected as methylhippuric acid in urine. The following tests are available:

Material: 10 ml urine

Preanalytics: cooled dispatch

TAT: 7-10 days*

Method: GCMS

Units: µg/L

Ref.- range: < 1500

For complete list of laboratory test offered at Freiburg Medical Laboratory, please visit <http://www.fml-dubai.com/parameter-listings/>